



- ① 调节螺栓 (PP)
- ② 活塞底座 (PP)
- ③ 外壳 (PP)
- ④ 容量设定螺栓 (PP, PEEK)
- ⑤ 针筒活塞(PFA、ETFE)
- ⑥ 活塞外壳 (PP)
- ⑦ 玻璃活塞壳 (Borosilicate glass)
- ⑧ 阀组件 (PFA), 外壳 (PP)
- ⑨ 密封球 (Borosilicate glass)
- ⑩ 进气阀 (ETFE/PTFE)
- ⑪ 伸缩式进液管 (FEP/ETFE)
- ⑫ 后盖 (PP)
- ⑬ 阀门开关 (PP)
- ⑭ 排液管套管 (PP)
- ⑮ 排液/循环阀(PFA/PTFE)
- ⑯ 排液管固定螺栓 (PP)
- ⑰ 排液管 (FEP)
- ⑱ 排液管盖 (PVDF)
- ⑲ 附属工具 (Size2.5-10ml, size25-100ml)
- ⑳ 止动环(PP)
- ㉑ 固定夹