

## Autoclavable Bottle-top Dispenser

# DISPET EX

- Autoclavable. (121°C for 20 minutes)
  - 3 sized adaptors for fitting with various types of bottles.
  - Chemical-guard materials used for liquid contact parts.
  - Adjustable Inlet tube for various types of bottles.
  - Liquid returning mechanism for minimizing the waste of reagents.
- \*Bottles are available as optional parts.



### Specifications

Cat. No	Volume range(mL)	Increment(mL)	Accuracy(%)	Precision(%)	Thread	Attached adapters
00-DPX-25	0.25 – 2.5	0.05	±0.6	≦0.1	GL-32	GL-28, GL-38, GL-45
00-DPX-50	0.5 – 5.0	0.1	±0.5	≦0.1	GL-32	GL-28, GL-38, GL-45
00-DPX-100	1.0 – 10.0	0.2	±0.5	≦0.1	GL-32	GL-28, GL-38, GL-45
00-DPX-250	2.5 – 25.0	0.5	±0.5	≦0.1	GL-45	GL-28, GL-32, GL-38
00-DPX-500	5.0 – 50.0	1	±0.5	≦0.1	GL-45	GL-28, GL-32, GL-38
00-DPX-1000	10.0 – 100.0	2	±0.5	≦0.1	GL-45	GL-28, GL-32, GL-38

## Bottle-top Dispenser

# DISPET Mini

- Suitable for repetitively micro volume dispensing.
  - Light and compact model for easy-to-carry-out.
  - Chemical-guard materials used for liquid contact parts.
- \*Bottles are available as optional parts.



### Specifications

Cat. No	Volume range (μL)	Accuracy (%)	Precision (%)	Thread
00-DPM-250	100 / 250	±2.0	≦0.4	GL-28
00-DPM-1000	500 / 1000	±1.0	≦0.2	GL-28
00-DPM-2000	1000 / 2000	±1.0	≦0.2	GL-28