

Fixed Volume and Triple Volume Micro Pipette

# Nichipet F&V

- The light and slender design along with the smooth movement of the plunger reduces the operator's fatigue.
- Built-in tip ejector.
- PTFE (Fluorine resin) is used in the airtight chamber of the instrument, so the instrument keeps its air tight seal and precise reproducibility for long hours.



**Nichipet F**

- ▷ Suitable when only one volume is needed.
- ▷ A wide range of standard volumes from 10µL – 1000µL.
- ▷ Other volumes within the range available for special order.

**Nichipet V**

- ▷ Three different volume models.
- ▷ Volume change by simply turning the knob.

**Specifications**

Type	Cat. No	Volume setting	Volume range(µL)	Accuracy(%)	Precision(%)	Usable tips
F	00-NPF-10	Fixed	10	±1.2	≤0.4	BMT-SS/SSR, CT-SS/SSR
	00-NPF-20		20	±0.8	≤0.3	
	00-NPF-25		25	±0.8	≤0.3	
	00-NPF-50		50	±0.6	≤0.3	BMT-SE/SER, SG/SGR, CT-SE/SER
	00-NPF-100		100	±0.6	≤0.3	
	00-NPF-200		200	±0.6	≤0.2	
	00-NPF-500		500	±0.5	≤0.2	
	00-NPF-1000		1000	±0.5	≤0.2	BMT-L/LR, CT-L/LR
V	00-NPV-S	Three variable settings	10/20/50	±1.2 – ±0.7	≤0.6 – ≤0.3	BMT-SE/SER, SG/SGR
	00-NPV-M		50/100/200	±0.7 – ±0.5	≤0.4 – ≤0.3	BMT-SE/SER, SG/SGR, CT-AG/AGR, SG/SGR
	00-NPV-L		200/500/1000	±0.7 – ±0.5	≤0.3 – ≤0.2	BMT-L/LR, CT-L/LR

\*The AC and CV values are the values by use of the disposable tips described in the catalogue.  
 \*See P15-16 for tip selection.

## Accessories

**Stand**



Cat. No ..... 00-MLT-STD2(ABS)  
 Storing capacity ..... 6pcs  
 Compatibility ..... Nichipet Premium  
 Nichipet EXII  
 Nichipet EX PlusII  
 Nichiryo Le



Cat. No ..... 00-BM-STD(Acrylic resin)  
 Storing capacity ..... 6pcs  
 Compatibility ..... Nichipet Premium  
 Nichipet EXII  
 Nichipet EX PlusII  
 Nichipet F&V



Cat. No ..... 00-NLE-STD(ABS)  
 Storing capacity ..... 1pcs  
 Compatibility ..... Nichipet Premium  
 Nichipet EXII  
 Nichipet EX PlusII  
 Nichiryo Le



Cat. No ..... 00-7000S(PVC)  
 Storing capacity ..... 2pcs  
 Compatibility ..... Nichipet 7000

**Tray**



Cat. No ..... 00-88-TY  
 Q'ty ..... 50pcs  
 Compatibility ..... Nichipet 7000

**Filter**

**For NPP1000**

Cat. No ... 00-NP-0601000  
 Q'ty ..... 10pcs

**For NPP5000**

Cat. No ... 00-NP-0605000  
 Q'ty ..... 100pcs

**For NPP10000**

Cat. No ... 00-NP-0610000  
 Q'ty ..... 100pcs

**Tube Opener**



Cat. No ..... 00-NTO  
 Q'ty ..... 3pcs  
 Compatibility ..... Nichipet Premium  
 Nichipet EX II  
 Nichipet EX Plus II

**For NPXII 1000**

Cat. No ... 00-NX2-0601000  
 Q'ty ..... 10pcs

**For NPLOII 1000**

Cat. No ... 00-LO2-0601000  
 Q'ty ..... 10pcs

**For 5000,10000**

Cat. No ..... 00-NX2-0615000  
 Q'ty ..... 100pcs  
 Compatibility ..... NPP, NPX2, NPLO2