


Disposable Tips

MICRO PIPETTE TIP

- The fine design and finishing of the top of the tips reduce the surface tension of liquid and ensure the tip to be residue free.
- The excellent sealing and discharge ability enable high accuracy of pipetting.
- Chemical and heat resistant, for which the advanced material "Highly purified polypropylene" is used.
- Fully autoclavable(121°C for 20 minutes). Racks can be autoclaved as well. (*Exceeding the recommended standard for autoclaving may cause deformities)
- The tips being able to be attached directly from the racks avoids dirt and dust from getting in, thus rising the efficiency of the operations.



Specifications

Type	Cat. No.	Specification	Volume range(μL)	Color	Length of tips	Q'ty(box)	Usable models	
UT	00-BMT-UT	Bulk tip	0.1 – 10	Clear	31.4mm	1000	NPP, NPX2, NPLO2-2, NP7-8V, 12V	
	00-BMT-UTRB	Racked tip				960(96pcs.x10 ⁶ /s)		
	00-BMT-UTRB	Racked tip (sterilized)				1000		
	00-CT-UT	Bulk filter tip				960(96pcs.x10 ⁶ /s)		
	00-CT-UTR	Bulk filter tip (sterilized)				1000		
SS	00-BMT-SS	Bulk tip	0.5 – 10	Clear	45.7mm	1000	NPP, NPX2, NPLO2, NLE, NPF-10	
	00-BMT-SSR	Racked tip				960(96pcs.x10 ⁶ /s)		
	00-BMT-SSRB	Racked tip (sterilized)				1000		
	00-CT-SS	Bulk filter tip				960(96pcs.x10 ⁶ /s)		
	00-CT-SSR	Bulk filter tip (sterilized)				1000		
SG	00-BMT-SG	Bulk tip	2 – 200	Clear	53mm	1000	NPP, NPX2, NPLO2, NLE-20, 100, 200, NP7-8S, 8L, 12S, 12L, NPV-S, M, NPF20 – 200	
	00-BMT-SGR	Racked tip				960(96pcs.x10 ⁶ /s)		
	00-BMT-SGRB	Racked tip (sterilized)	2 – 20			1000		
	00-CT-SG	Bulk filter tip				960(96pcs.x10 ⁶ /s)		
	00-CT-SGR	Bulk filter tip (sterilized)				1000		
SE	00-BMT-SE	Bulk tip	2 – 200	Clear	51.1mm	1000	NPX2, NPLO2, NLE-20, 100, 200, NPV-S, M, NPF20 – 200	
	00-BMT-SER	Racked tip				960(96pcs.x10 ⁶ /s)		
	00-BMT-SERB	Racked tip (sterilized)	10 – 100			1000		
	00-CT-SE	Bulk filter tip				960(96pcs.x10 ⁶ /s)		
	00-CT-SER	Bulk filter tip (sterilized)				1000		
AG	00-CT-AG	Bulk filter tip	20 – 200	Clear	66mm	1000	NPP, NPX2, NPLO2, NLE-200, NPV-M	
	00-CT-AGR	Bulk filter tip (sterilized)				960(96pcs.x10 ⁶ /s)		
	00-CT-AGR	Bulk filter tip (sterilized)				1000		
K	00-BMT-K	Bulk tip	50 – 300		Clear	58.9mm	1000	NP7-8K, 12K
	00-BMT-KR	Racked tip					960(96pcs.x10 ⁶ /s)	
	00-BMT-KRB	Racked tip (sterilized)		1000				
L	00-BMT-L	Bulk tip	100 – 1000	Clear		77mm	1000	NPP, NPX2, NPLO2, NLE-1000, NPV-L, NPF500, 1000
	00-BMT-LR	Racked tip					1000(100pcs.x10 ⁶ /s)	
	00-BMT-LRB	Racked tip (sterilized)			1000			
	00-CT-L	Bulk filter tip			1000(100pcs.x10 ⁶ /s)			
	00-CT-LR	Bulk filter tip (sterilized)			200			
X	00-BMT-X	Bulk tip	1000 – 5000	Clear	132mm	100(50pcs.x2 ⁶ /s)	NPP, NPX2, NPLO2-5000, 11DG-5000	
	00-BMT-XR	Racked tip				200		
	00-BMT-XRB	Racked tip (sterilized)				100(50pcs.x2 ⁶ /s)		
Z	00-BMT-Z	Bulk tip	1000 – 10000		Clear	158mm	200	NPP, NPX2, NPLO2-10000

* The tips are available for pipette models with in the volume range of the tip mentivnal above.