

Glass Tip Pipette

Nichipet ECO

- The accuracy and precision of a micropipette with the stability and reliability of glass tips.
- The glass tip enables bias free liquid handling for endocrine disruptor's analysis.
- Four different sizes covering a wide volume range. (20μL to 10mL)
- Glass tip is secured to the unit with a double O-ring design.
- Comes with a PTFE nozzle and ceramic piston.



Glass Tip




Specifications

| Cat. No | Volume range(μL) | Increment(μL) | Accuracy(%) | Precision(%) | Usable tips |
|-------------|------------------|---------------|-------------|--------------|-------------|
| 00-NPC-200 | 20 - 200 | 0.2 | ±6.0 - ±3.0 | ≤6.0 - ≤3.0 | Glass Tip S |
| 00-NPC-1000 | 100 - 1000 | 1.0 | ±5.0 - ±3.0 | ≤5.0 - ≤3.0 | Glass Tip L |
| 00-NPC-5000 | 1000 - 5000 | 10.0 | ±5.0 - ±3.0 | ≤5.0 - ≤3.0 | Glass Tip X |
| 00-NPC-10ML | 1000 - 10000 | 50.0 | ±6.0 - ±3.0 | ≤6.0 - ≤3.0 | Glass Tip Z |

● Glass tips

| Cat. No | Volume range(μL) | Q'ty(box) |
|----------|------------------|-----------|
| 00-GLT-S | 20 - 200 | 10 |
| 00-GLT-L | 100 - 1000 | 10 |
| 00-GLT-X | 1000 - 5000 | 10 |
| 00-GLT-Z | 1000 - 10000 | 10 |

Further research into "Environmental Hormone" is strongly requested to clarify the precise cause and effect on the human system and wild living things. Glass tips rather than standard plastic pipette tips are more suitable for this research.



Variable Volume Macro Pipette

MODEL 3100

- Four different volume models from 1 mL to 10mL.
- Light weight and slender body.
- Easy to operate.
- Liquid touches only disposable tip.
- PTFE piston ensures complete sealing of the syringe.



Specifications

| Cat. No | Volume range(mL) | Increment(mL) |
|-----------|------------------|---------------|
| 00-31N-1 | 0.2 - 1.0 | 0.1 |
| 00-31N-2 | 0.4 - 2.0 | 0.1 |
| 00-31N-5 | 1.0 - 5.0 | 0.5 |
| 00-31N-10 | 2.0 - 10.0 | 0.5 |

- Accuracy and Precision ± 1.0% / ≤ 1.0%
- Disposable tips

| Cat. No | Volume range(mL) | Cloror | Q'ty(bag) |
|----------|------------------|--------|-----------|
| 00-MPT-1 | 0.2 - 1.0 | Clear | 1000 |
| 00-MPT-2 | 0.4 - 2.0 | Yellow | 200 |
| 00-BMT-X | 1.0 - 5.0 | Green | 200 |
| 00-BMT-Z | 2.0 - 10.0 | Clear | 200 |

Variable Volume Ultra Macro Pipette

MODEL 5100

- Two models, each with three standard volume settings.
- Volume change done by placing a volume spacer in the pipette.
- Any special volume spacer between 15 to 25 mL for 5100S and 30 to 50 mL for 5100L is available for special order.



Specifications

| Cat. No | Volume setting(mL) | | |
|----------|--------------------|------|------|
| 00-5100S | 15.0 | 20.0 | 25.0 |
| 00-5100L | 30.0 | 40.0 | 50.0 |

- Accuracy and Precision ± 0.3% / ≤ 0.3%
- Disposable tips

| Cat. No | Volume range(mL) | Q'ty(bag) |
|---------|------------------|-----------|
| 00-51TS | 25.0 | 5 |
| 00-51TL | 50.0 | 5 |